

Vesicles

Vesicular transport assay

Vesicular transport assay

With the Cell4Pharma Vesicular Transport Assay substrate determination, inhibition studies and several kinetic studies can be performed. Transport can be measured radioactive, fluorescent or analytical (LCMS).

Equipment:

- Plate Shaker
- Water bath
- Multiscreen HTS-Vacuum Manifold filtration device (Merck Millipore)
- Multichannel pipets (advised) Scintillation counter, fluorescent detector or LC-MS/MS
- Scintillation counter, fluorescent detector or LC-MS/MS

Disposables:

- SV-shape 96-wells incubation plates (651101, Greiner-bio-one) (make hole in corner to remove air and get better water contact)
- Filter plates
 - PVDF (MSHVN4B50, Merck Millipore)
 - Glass Fiber (MSFBN6B50, Merck Millipore)

Radioactive measurement:

- Punching tips (MADP19650, Merck Millipore)
- Counting vials
- Topcount adapter, **for 96-wells counting only** (MSTPCWH50, Merck Millipore)

Fluorescent measurement:

- Black flat bottom 96-wells plate (655096, Greiner-bio-one)

Analytical measurement:

- Flat bottom 96-wells plate (655101, Greiner-bio-one)

Buffers/Solutions:

- Scintillation fluid (6013199, Perkin Elmer)
- Cell4Pharma Stop/Wash buffer (2x)
- Cell4Pharma assay buffer (2x): 500 mM Sucrose, 20 mM Tris-HEPES, pH 7.4
- MgCl₂, 500 mM
- AMP, 100 mM pH 7.4
- ATP, 100 mM pH 7.4

Vesicular transport assay

Assay suggestions:

- Performing incubations in triplicate
- For every test condition, dedicate three wells for AMP-containing buffer (background signal) and three wells for ATP-containing buffer. The difference between AMP- and ATP-containing conditions represents ATP-dependent uptake

Preparation:

- Prepare substrate mixture according to the following scheme per well:

Mixture for each well:		
Cell4Pharma assay buffer 2x	12.5	µl
Cell4Pharma MgCl ₂ (500 mM)	0.6	µl
Substrate	x *	µl
Inhibitor	x *	µl
Cell4Pharma AMP or ATP (100 mM)	1.2	µl
MilliQ	x *	µl
		Total 25 µl

Maximum final concentration of DMSO 1%

Procedure:

- Place the v-shaped 96-wells plate on ice.
- Thaw the vesicles at 37°C for 1 minute and place on ice.
- Pipet 5 µl vesicles at the bottom of the well.
- Add 25 µl/well ice cold substrate mixture prepared earlier (mix on shaker for 10 sec).
- Incubate the plate in a 37°C water bath for x minutes (transporter dependent)
- Dilute the original Cell4Pharma Stop/Wash buffer 2X first and stop the reaction by transferring the plate on ice and immediately add 150 µl ice-cold diluted Cell4Pharma Stop/Wash buffer
- Prewet the 96-wells filter plate # with 200 µl Cell4Pharma Stop/Wash buffer and apply vacuum to remove buffer just before transferring the samples.
- Transfer the samples to the filter plate and wash twice with 200 µl Cell4Pharma Stop/Wash buffer by applying vacuum. Be sure the filters are dry after the second wash step.
- **The following step is dependent on the method of sample analysis:**
 - **a.** Radioactivity counting in vials: remove plastic bottom underneath the filter plate and transfer filters using punching tips and add scintillation liquid
 - **b.** 96-wells radioactivity counting: carefully remove plastic bottom underneath the filter plate, place a topcount adapter and add scintillation liquid
 - **c.** Analytical measurements: place a flat-bottom 96-wells plate in the vacuum device, extract the samples with solvent (substrate dependant) by applying vacuum.
 - **d.** Fluorescent measurements: place a black flat-bottom 96-wells plate in the vacuum device, extract the samples with solvent (substrate dependant) by applying vacuum.